

deltaPREP Blood DNA Mini Kit (MDX)

The deltaPREP Blood DNA Mini Kit (MDX) is a highly efficient tool for directly isolating DNA from whole blood samples of up to 400 µl.

The optimized kit for extraction of DNA from blood guarantees extremely high-quality DNA that can be used immediately for photometric determinations, PCR or other downstream applications. The deltaPREP Blood DNA Mini Kit (MDX) is very fast and easy to use. Isolation of genomic DNA – free of inhibitors or impurities – can be completed in just 24 minutes.

Product Name: deltaPREP Blood DNA Mini Kit (MDX)

Product details

Extract: DNA

Reactions: 10 / 50 / 100 / 250 (either or)

Sample type/Starting material: Blood, fresh or frozen (RNA only fresh)

Specifications:
Starting material

- Whole blood samples (up to 400 µl)
- Fresh or frozen blood
- Stabilizers: EDTA or citrate

Average yield

Depends on sample and the used volume
Up to 30 µg gDNA

Extraction time

Approx. 24 minutes, including lysis

Binding capacity

Column binding capacity: > 50 µg gDNA

Average purity

1,7 - 2,0

The online shop

Price: € 42.00

Content: 10 reactions

Please select packing

10 reactions ▼