

# innuSPEED Lysis Tube J

**Lysis Tube J** is a 2.0 mL tube with screw cap (blue) containing optimized steel beads (4.7 mm) for efficient disruption of hard samples as wood, seed or rice grain using SpeedMill or other homogenizers.

## Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 or 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

**Product Name:** innuSPEED Lysis Tube J

## Product details

**Lysis Tubes:** Steel beads 4.7 mm

## The online shop

Price: € 31.50

Content: 10 reactions

Please select packing

10 reactions ▼